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Thank you for choosing STERN PODIA products.

Before operating the unit, we invite you to read this manual carefully and keep it in a safe place: this will enable you to discover and optimally exploit all the potentialities of your equipment and ensure that constant levels of performance, safety and reliability are maintained over time.

And, whatever circumstances may arise, you can count on the professionalism of your S.M.D. dealer, the expertise of its technicians and our reliability as constructors.

**We wish you all the best for your work!**

<b>• TABLE OF CONTENTS</b>	
TECHNICAL CHARACTERISTICS	Pg.9
WARNINGS AND CAUTIONARY NOTES	Pg.10
INSTALLATION AND CONNECTIONS	Pg.11
OPERATION	Pg.11
PROGRAMS AND PROGRAMMING	Pg.12
ROUTINE MAINTENANCE	Pg.14
SPECIAL MAINTENANCE	Pg.14
CLEANING AND DISINFECTION	Pg.14
LIST OF ACCESSORIES INCLUDED	Pg.14
CONFORMITY DECLARATION	Pg.17

### T-FORM TECHNICAL CHARACTERISTICS

Power supply	230 VAC single phase /50Hz
Absorption	1700 VA
Allowed mains voltage fluctuation	+/- 10%
Mains fuses	n.2 (T.3,15A/250V)
Liquids protection rating	IPX 0
Operating temperature range	10-35° C
Storage and transport temperature	-10°/70°C
Operating humidity range	45/90 %
Storage and transport humidity	10/95 %
Working atmospheric pressure	500 / 1060 mbar
Noise level	50 dB (A)
Thermoformer weight	15 Kg
Vacuum pump weight	4 Kg
Work surface dimensions	285 mm x 340 mm
Thermoformer dimensions	520 mm x 520 mm x130 mm (h)
Vacuum pump dimensions	194 mm x 320 mm x 182 mm (h)

## WARNINGS AND CAUTIONARY NOTES

### Intended use:

The T-FORM thermoformer is intended to be installed in podologic cabinets and orthopaedic centres and may be used solely by qualified personnel for the thermoforming of arch supports using a suitable material.

Do not remove the labels affixed to the device. Should they deteriorate, have them replaced.

Compliance with directive 2006/42/CEE, 2014/35/UE e 2014/30/UE and subsequent amendments.

### Preparations to be made by the user:

The user must arrange for installation in a suitable room equipped with an electrical system complying with current regulations.

The device should be installed in a dry, well-lit environment complying with the requirements of current laws.

By current laws/regulations we mean the regulatory framework in force in the user's country.

Make sure that the power supply line of the units is fitted with an external thermal magnetic differential circuit breaker with a rating of 16A and trip threshold of 10 mA.

### Climate:

In extreme climate conditions (heat, cold, humidity), after removing the unit from the packing container it is advisable to wait a few hours before switching it on: this precaution will allow time for the evaporation of any condensate that may have formed inside the packing container.

### Recommendations for safeguarding the environment

The different packing materials must not be discarded in the environment but must rather be disposed of in accordance with the regulations established by the local authorities.

Scrapping must be carried out according to local waste disposal regulations. Contact your local Stern Podia dealer.

Limit the operation of the heating surface to the time that is strictly necessary to complete the planned operations.

### General warnings and cautionary notes

Any repairs or servicing must be performed solely by a qualified technician.

Avoid exposing any part of the unit to atmospheric agents.

Do not place weights or sharp or pointed objects on the unit, as this may cause damage to the silicone membrane.

Keep the part in contact with the silicone membrane clean in order to prevent deterioration of the membrane itself and guarantee correct sealing of the thermoforming compartment.

Do not pour water directly onto the unit, either when open or closed, for any reason. To facilitate cooling, wipe a moist cloth over the outer part of the silicone membrane.

Avoid contact between the unit and wet parts of the body and do not use the unit while barefooted.

Do not place flammable materials inside the unit.

Make sure that power cables or accessory parts of the unit do not come into contact with the heated parts of the unit itself.

Daily, check the water filter fit on the vacuum pump.

### Precautions

The grinding box needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information contained herein.

The portable devices may affect the operation of the device.

Use only accessories that are specified in this document to ensure the maintenance of the EMC performance of the device.

Do not use together with other nearby devices to avoid mutual interference.

Attention to the connection of the electrical power strip leads to the creation of a SYSTEM EM, since this situation may result in a reduced level of security. About the applicable requirements SYSTEM EM the ORGANIZATION

RESPONSIBLE must refer to this standard.

Upon request S.M.D. provides circuit diagrams, the parts list, descriptions, calibration instructions or other information to assist the technical support personnel during installation or repair.

Please note if you find any fissures or cracks on the power cord or the same was partly cut off, operate the switch activity by service personnel qualified by SMD

### **WEAR SUITABLE SAFETY GLOVES**

Do not use materials for thermoforming if their intended use is unknown.

**S.M.D. cannot be held liable for any damages deriving from failure to observe the warnings and instructions contained in this handbook.**


## INSTALLATION AND CONNECTIONS

Connect the air hose to the connector (2) provided on the rear. Connect the power supply cable of the thermoformer (3) to the mains power supply. The green lighted power on/off switch (4) on the rear of the unit and the control panel display turn on to indicate that the device is ready for operation.

## OPERATION

The T-Form thermoformer represents the state of the art in the podiatric practice of arch support moulding. Thanks to the electronic microprocessor control, the operator can adjust different parameters based on his/her specific needs and according to the different types of materials and relevant thicknesses. Each chosen setting is shown on the control panel display. Five customisable automated programs are available along with two more programs, also customisable, dedicated to specific impression-taking techniques.

+ and – buttons for adjusting the temperature of the work surface from 40°C to 240°C

+ and – buttons for adjusting the vacuum level from 0 to -998 mbar

+ and – buttons for setting the time that the set temperature (in the TERMO mode) and vacuum level (VACUUM mode) will be maintained



Button for selecting the desired program. If it is kept pressed down for around 3 seconds, the program parameter settings shown on the display will be memorised

Button for starting the selected program.	
Button for stopping the selected program or interrupting a cycle in progress. Sequential button for switching from the TERMO MODE to the VACUUM MODE and vice versa within each program. (excluding program A)	

### PROGRAMS AND PROGRAMMING

<p><b>PROGRAMS FROM 1 TO 5</b></p> <p>By pressing the PRG button repeatedly in sequence you can select one of the programs from <b>1 to 5 shown on the display</b>.</p>	
<p>For each of these programs it is possible to set;</p> <p>in the TERMO mode</p> <ul style="list-style-type: none"> <li>- thermoforming temperature;</li> <li>- time for which the thermoforming temperature will be maintained;</li> </ul> <p>and, by switching over to the VACUUM mode</p> <ul style="list-style-type: none"> <li>- vacuum level to be used for thermoforming;</li> <li>- time for which the set vacuum level will be maintained once it has been reached;</li> </ul>	
<p>If you keep the PRG button pressed down for approximately 3 seconds, an audible signal will be emitted to confirm memorization of the parameters selected for the program indicated on the display.</p> <p>From this moment on the unit will maintain the programming even after being turned off and back on again by means of the main switch.</p> <p>In any case, you can edit the parameters at any time by means of the appropriate buttons. Any changes will be maintained, if not memorised, until you switch over to a different program.</p>	
<p>Start the thermoforming cycle by pressing the START button. The unit will heat up the work surface and on reaching the set temperature will begin the countdown of the time for which the set temperature will be maintained. An audible signal and the message "OK" flashing on the display will indicate to the operator when this time has elapsed. In this status the unit maintains both the vacuum and the temperature modes.</p> <p>Interrupt the TERMO mode by pressing the STOP button to enable the lid to open; the unit will automatically switch into the VACUUM mode.</p> <p><u>Warning! If necessary, press the STOP button again to repeat the TERMO cycle.</u></p> <p>Activate the VACUUM mode by pressing the START button. Upon reaching the set vacuum level the unit will begin the countdown until the interruption of the cycle, which will enable the lid to open. The unit will automatically switch into the TERMO mode.</p> <p><u>Warning! If necessary, press the STOP button again to repeat the VACUUM cycle.</u></p>	

**PROGRAM "A"**

This is a dedicated program for taking impressions directly using the accessories provided.

Connect the hose provided to one of the quick-release connectors on the front of the unit and the vacuum level remote control to the socket provided.

Set the parameter for the minimum vacuum level; in this manner, when changes are made via the remote control it will not be possible in any case to go below the set value.

By keeping the PRG button pressed down for approximately 3 seconds you can memorise the parameters selected for the program indicated on the display.

Warning! In this program the TERMO mode is disabled.

**PROGRAM "B"**
**PROGRAMMING**

Select:

TEMPERATURE: for example, 110° (depending on material type)

WARMING TIME: this operation is only necessary when the plate is cold.

Press  to select the « VACUUM » mode:

MAXIMUM VACUUM: 100% (not changeable)

MINIMUM VACUUM: adjustable up to 5%

To store the set data, press  until you can hear a confirmation BIP

**USE**

1-Press  to exhaust the pads

2-The patient can get up on the pads made rigid by the vacuum.

3-Open the pads and let air go in by pressing at the same time the RED and GREY buttons on the remote control until you can hear a confirmation BIP.

4-Taking the footprint: exhaust the pads by pressing the GREY button until you can hear a confirmation BIP.



5-Deaccrease the vacuum down to the minimum set to correct the footprint: press the RED button until you can hear a confirmation BIP.

6-Correct the footprint.

7-After correcting the footprint, press the GREY button to confirm the corrections

8-The patient can now get off the pads.

9-Thermo-forming of insoles under patient's feet.

**ROUTINE MAINTENANCE**

After each thermoforming cycle, thoroughly clean the work surface and silicone membrane with a damp cloth, so as remove any scraps or traces of glues.  
 Periodically check the oil level in the vacuum pump.  
 Check the silicone membrane of the lid and protective guard of the work surface for wear. If necessary, request a replacement from a service centre.  
 Daily, check the water filter fit on the vacuum pump.

**SPECIAL MAINTENANCE**

For any special maintenance work, contact an authorised service centre.  
 Precisely describe the fault or defect detected to the service engineers; also specify the serial number of the bench, shown on the label affixed on the rear of the unit.

**CLEANING AND DISINFECTION**

For the cleaning and disinfection of exterior surfaces use standard commercially available disinfectants, checking that the active ingredients listed are present in the following maximum concentrations:

- Glutaraldehyde 2%:10%
- Ethanol 96%:40%
- Formaldehyde:0.01%
- Glyoxal:0.15%
- Propanol:35%

**WARNING!**

Do not exceed the specified concentrations: do not use products containing alcohol, ammonia, abrasive substances or benzene; wipe a cloth dipped in water over disinfected parts (disinfectants may damage surfaces even if diluted).  
 The products should be applied using a soft cloth dipped in disinfectant: to avoid infiltrations that may impair performance, do not spray the disinfectant directly on the unit.  
 S.M.D. disclaims all liability for any injury or damage resulting from inobservance of the warnings and instructions contained in this manual.

**LIST OF ACCESSORIES INCLUDED**

<b>1</b>	HOSE FOR CONNECTION TO BAG FOR DIRECT IMPRESSION
<b>2</b>	TIE FOR CLOSING BAG DURING VACUUM CYCLE
<b>3</b>	FITTING FOR CONNECTION TO SILICONE BAGS
<b>4</b>	BAGS FOR DIRECT IMPRESSION OF FOOT
<b>5</b>	REPLACEMENT FUSES (fuse holder situated on rear of T-FORM)
<b>6</b>	REMOTE CONTROL FOR VACUUM

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