---

**D**

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---

**GB**

Thank you for choosing STERN PODIA products.

The equipment you have purchased has been designed and manufactured to meet your present needs and, thanks to its in-built versatility and adaptability, your future needs too.

Before using this unit, it is good practice to read this handbook carefully and keep it safely.

Doing so will provide you with a full understanding of the unit, allow you to make the best use of its potentiality and maintain top-level performance, safety, and reliability throughout its working life.

And, whatever circumstances may arise, you can count on the professionalism of your S.M.D. dealer, the expertise of its technicians and our reliability as constructors.

**We wish you all the best for your work!**

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### GB TABLE OF CONTENTS

TECHNICAL SPECIFICATIONS	Page 7
KEY TO SYMBOLS	Page 8
WARNINGS	Page 12
CLEANING AND DISINFECTION	Page 13
IDENTIFICATION - DIMENSIONS	Page 14
GENERAL INSTRUCTIONS	Page 15
CONTROL PANEL	Page 16
CHAIR CONTROL	Page 17
INSTRUMENTS	Page 18
SELECTION OF INSTRUMENTS MAXIMUM ROTATION SPEED	Page 22
SPRAY FUNCTION	Page 23
FOOT CONTROL	Page 24
ROUTINE MAINTENANCE	Page 25
TROUBLESHOOTING	Page 27
SPARE PARTS LIST	Page 27
CONFORMITY DECLARATION	Page 28
WARNINGS ON ELECTROMAGNETIC INTERFERENCES	Page 31

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
**GB TECHNICAL SPECIFICATIONS**

Power supply	230 V AC single phase/50Hz
Podiatry unit absorption	70 VA max
Allowed mains voltage fluctuation	+/- 10%
Mains fuses	n.2 (T.6,3 A /250V)
Electrical protection rating	Class I -Type B
Liquids protection rating	IPX 0
Surgery electrical system	IEC 60364
Separate on-off switch	Capacity 16A /trip threshold 10A
Air inlet	In feed pressure 5.0 – 7.0 bar Flow-rate 10 l/min with Turbine 30 l/min Connection pipe diam.6x8
Suction	Suction devices powered directly from the unit: max limit of absorption, 2000 VA
Spray liquid tank capacity	1 litre
Podiatry unit top max load capacity	Kg.15
Operating temperature range	10 / 35 °C
Storage and transport temperature	-10 / 70 °C
Operating humidity range	45 / 95 %
Storage and transport humidity	10 / 90 %
Atmospheric pressure	860 / 1060 mbar
Unit weight	Kg.80 (with drawer's element)
Unit dimensions (cm)	790 x 640 x 920(H) mm
Drawer unit load capacity	Kg.3,00

Weights are indicative only and may vary according to the exact configuration.

<b>D</b> <b>VERWENDETE</b> <b>SYMBOLS</b>	<b>E</b> <b>SÍMBOLOS</b>	<b>GB</b> <b>KEY TO SYMBOLS</b>	
Wechselstrom	Corriente alterna	Alternating Current	
Schutzerdung	Conexión a tierra	Earth connection	
Gerät Klasse I - Typ B	Aparato de clase 1 Tipo B.	Class 1 Type B Equipment	
Hochspannung (230 V)	Alta tensión (230 V)	High Voltage (230 V)	
Eingeschaltet (an das Stromnetz angeschlossen)	Activado (conectado a la red de alimentación)	On (Connected to power supply)	
Ausgeschaltet (an das Stromnetz nicht angeschlossen)	Desactivado (desconectado de la red de alimentación)	Off (disconnected from power supply)	
Lesen Sie die Gebrauchsanweisung	Leer el manual del usuario	See the operative instructions	
ACHTUNG, vor Gebrauch das Handbuch lesen	ATENCIÓN, consultar el manual de instrucciones antes de la utilización.	WARNING: see instruction manual before use	
Gemäß der Richtlinie 2007/47/EEC "MEDICAL DEVICES" "Medizinprodukte"	En conformidad con la directiva 2007/47/EEC "MEDICAL DEVICES"	Conforms to EEC Directive 2007/47/EEC "MEDICAL DEVICES"	
Technischer Kundendienst (nur von Fachpersonal ausgeführt)	El Servicio postventa (intervención realizada sólo por personal técnico calificado)	Customer Support Service (task to be done by qualified technical staff only)	<b>S.P.V.</b>
Im Autoklav sterilisierbar	Parte que se puede esterilizar en autoclave	Piece than can be sterilized in autoclave	

Entsorgen gemäß den anzuwendenden Gesetzen. Wenn Sie Fragen oder Bedenken haben, kontaktieren Sie bitte SMD	Eliminar de acuerdo con la ley aplicable. Si usted tiene alguna pregunta o duda por favor póngase en contacto con SMD	Dispose of in accordance with applicable law. If you have any questions or concerns, please contact SMD	
Sitzen verboten	Prohibido sentarse	No sitting	
Steigen verboten	Prohibido subir	Do not stand	
Diese Seite nach oben	Lado superior	Upper side	↑↑
Zerbrechlich	Frágil	Fragile	
Vor Nässe schützen	Proteger de la lluvia	Protect from rain	

Die auf dem Gerät erscheinenden Symbole entsprechen der IEC-Norm 601.1 und der ISO-Norm 9687.  
 Los símbolos que aparecen en el equipo cumplen con las normas IEC 601.1 e ISO 9687.  
 The symbols that appear on the equipment comply with IEC 601.1 and ISO 9687 standards.

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---

**GB WARNINGS****Purpose of the unit.**

The podiatry unit is intended for professional podiatrists as an aid to in-surgery operations that concern ailments of the foot and their treatment.

These units have been designed and constructed for use with patients of all ages.

Never the podiatry unit may be used in the presence of flammable anesthetic gas mixtures. Stern Podia units are not sterile devices.

The compartment with the "U.V." germicide lamp is not a sterilizing device.

When treating patients with pacemakers or similar devices the operator must take the necessary precautions: refer to the specific literature on the matter.

Do not exceed the loads indicated in the "technical characteristics" section. Do not place any additional loads on or plug any device to the instrument console or to the lamp.

Do not remove the labels on any unit/chair devices. Should these deteriorate replace them.

The device is not intended to be used in emergency rooms or operating theaters.

**Power supply**

Make sure that the unit/chair electrical power supply is powered via an external differential cut-out switch with a 16A current-carrying capacity and a 10mA trip threshold.

Make sure that unit power supplies are separate.

Make sure that the podiatry unit power supply is fitted with an efficient grounding (earth) system and that wiring complies with the standards in force.

Make sure that the compressed air feed line has an external close-off valve.

Make sure that the air hoses are bled carefully and thoroughly before connecting to the podiatry unit.

**Climate**

Under extreme conditions (heat, cold, humidity) it is advisable to let a few hours pass between unpacking the equipment and using it for the first time. This precaution provides the time needed to eliminate any condensation that may have formed inside the packaging.

**Precautions**

The unit needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information contained herein.

The portable devices may affect the operation of the device.

Use only accessories that are specified in this document in order to ensure the maintenance of the EMC performance of the device.

Do not use together with other nearby devices to avoid mutual interference.

Attention: the connection of the electrical power strip leads to the creation of a SYSTEM EM, and this situation may result in a reduced level of security. About the requirements applicable to a EM SYSTEM, the ORGANIZATION RESPONSIBLE must refer to this standard.

The unit must be connected to an external air supply, whose characteristics are given in the Technical Specifications section. For the insertion of the tube (according to DIN 74324 specifies 8X1 PA11 PHL) a quick connector of 6x8 is provided.

These operations must be carried out by technicians with appropriate knowledge of the device.

For the extraction, there are no special requirements.

To preserve the condition of the buffer battery present in the unit, make sure that it does not remain disconnected from the power supply for longer than three months.

The battery replacement should be carried out by technicians with appropriate knowledge of the device.

Upon request S.M.D. provides circuit diagrams, spare parts list, descriptions, calibration instructions or other information to assist the technical support personnel during installation or repair.

Should you find any fissures or cracks on the power cord or should it be partly cut off, get SMD qualified service personnel to change it.

**S.M.D. cannot be held liable for any damages deriving from failure to observe the warnings and instructions contained in this handbook.**

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<b>GB</b> <b>TROUBLESHOOTING</b>	
Turning main switch, does not turn on the unit.	<ul style="list-style-type: none"> <li>• Check the external differential cut-out switch.</li> <li>• Check the fuses.</li> <li>• Contact the CSS (Customer Support Service)</li> </ul>
Instruments are not enabled.	<ul style="list-style-type: none"> <li>• Check the main switch.</li> <li>• Contact the CSS.</li> </ul>
No water and/or air.	<ul style="list-style-type: none"> <li>• Check that pneumatic air selector is open.</li> <li>• Check that external feed lines are open.</li> <li>• Contact the CSS.</li> </ul>
Progressive reduction of fluid flow-rate.	<ul style="list-style-type: none"> <li>• Clean or replace the unit in feed filter cartridges (CSS)</li> </ul>
Progressive reduction of suction module/device performance.	<ul style="list-style-type: none"> <li>• Check the suction filter and replace if necessary.</li> <li>• Check cleaning the suction pipes and replace if necessary (CSS).</li> </ul>
Spray water is not fully nebulised.	<ul style="list-style-type: none"> <li>• Check that instrument spray nozzles are unobstructed; make any necessary adjustment.</li> <li>• Check water level inside the tank.</li> <li>• Contact the CSS</li> </ul>
Optic fibres do not function.	<ul style="list-style-type: none"> <li>• Check that the function has been enabled.</li> <li>• Contact the CSS</li> </ul>

	<b>ERSATZTEILLISTE</b>	<b>LISTA DE REPUESTOS</b>	<b>SPARE-PARTS LIST</b>
469600007	MIKROMOTOR MCX	MICROMOTOR MCX	MCX MICROMOTOR
469600006	MIKROMOTORSCHLAUCH	MANGUERA MICROMOTEUR	MICROMOTOR TUBE
478990005	TURBINESCHLAUCH	MANGUERA DE LA TURBINE	TURBINE TUBE
481042000	COMPLETE SYRINGE	JERINGA	SYRINGE
460233640	SCHUBLADENGRIF	TIRADOR DEL CAJON	DRAWER HANDLE
457604011	RAD MIT BREMSE	RUEDA CON FRENO	WHEEL WITH BRAKE

**DICHIARAZIONE DI CONFORMITA'**  
**DECLARATION DE CONFORMITE**  
**CONFORMITY DECLARATION**  
**KONFORMITATSERKLÄRUNG**

S.M.D. garantisce e dichiara sotto la propria responsabilità che i dispositivi medici /  
 S.M.D. déclare que les dispositifs médicaux  
 S.M.D. declares that the medical devices / S.M.D. erklärt, dass die Medizinprodukte

Tipo / Type	Marchio/Marque/ Trademark / Brandmarke	Modello / Modèle / Model / Modell	Matricola / <u>No. de</u> fabrication /Serial No. / Seriennummer. da/de/fm/von ÷ a/a/to/bis
Riunito podologico	Stern Podia	ATLAS	
Unité podologique			39103-00-0000 ÷ 39103-99-9999
Podiatric Unit			
Behandlungseinheit			

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Are medical devices of Class IIa, classification rule 9 (Encl. IX of the 93/42/EEC Directive modified by the 2007/47/CE).  
 They follow the fundamental safety requirements as indicated in the Encl. I of the 93/42/EEC Directive (as modified by the 2007/47/CE). The quality system implemented according to the Encl. II of the 93/42/EEC Directive is certified and evaluated by the Notified Body n.0476.

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## WARNINGS ON ELECTROMAGNETIC INTERFERENCES

The device has undergone EMC (Electromagnetic Compatibility) verifications for medical devices as per EN60601-1-2.

The device needs special precautions about the EMC and must therefore be installed and commissioned following the EMC information included in this manual.

Any portable and mobile radio communications equipment (including also transmission peripherals such as cables and antennas) could influence the functioning of the device. Therefore, they should be placed at a distance of at least 30 cm. Failure to observe this indication could result in a malfunctioning.

Furthermore, interferences could occur in proximity of any equipment marked with the following symbol

The use of cables different from the provided ones could cause an increase of the EMISSIONS or a decrease of the IMMUNITY of the device.

The device should not be used near other electrical or electronic equipment. If necessary, verifying that the device functions correctly. In case of faults or transient noises, not only the device stops working but it also clearly indicates that. The device can be used in any building, including also residential buildings and those directly connected to the public low-voltage power supply network which supplies buildings intended for residential use.

Based on the verifications carried out, the following further warnings are provided:

The device is suited for the use in a specific electromagnetic environment. The buyer or the end-user of the device should make sure that it is used in an electromagnetic environment having the following characteristics:

Emissions test	Conformity
Radiated or conducted RF CISPR 11	Group 1 Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations and flicker emission IEC 61000-3-3	Conform

Immunity test	Test level (and conformity test) EN60601-1-2 (maximum values)
ESD EN61000-4-2	± 8kV contact discharge ± 15kV air discharge
Electrical Fast Transient (EFT) EN61000-4-4	± 2kV power line ± 1kV I/O line
Surge EN61000-4-5	± 1kV differential mode voltage ± 2kV common mode voltage
Voltage dips, short interruptions and voltage variations EN61000-4-11	0% U <sub>n</sub> during 0.5 cycle 70% U <sub>n</sub> during 25 cycles 0% U <sub>n</sub> during 5 seconds
Power-frequency magnetic field (50/60 Hz) EN61000-4-8	30A/m
Radiated RF EN61000-4-3	3V/m 80MHz ÷ 2.7GHz
Conducted RF EN61000-4-6	3V 150kHz ÷ 80MHz