**I**

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**Bon travail !**

**GB**

Thank you for choosing STERN PODIA products.

The equipment you have purchased has been designed and manufactured to meet your present needs and, thanks to its in-built versatility and adaptability, your future needs too.

Before using this unit, it is good practice to read this handbook carefully and keep it safely.

Doing so will provide you with a full understanding of the unit, allow you to make the best use of its potentiality and maintain top-level performance, safety, and reliability throughout its working life.

And, whatever circumstances may arise, you can count on the professionalism of your S.M.D. dealer, the expertise of his technicians and our reliability as constructors.

**We wish you all the best for your work!**

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**GB TECHNICAL SPECIFICATIONS**

Power supply	230 V AC single phase/50Hz
Podiatry unit absorption	70 VA max
Permissible mains voltage fluctuation allowed	+/- 10%
Mains fuses	n.2 (T.6,3 A /250V)
Electrical protection rating	Class I -Type B
Liquids protection rating	IPX 0
Surgery electrical system	IEC 60364
Separate on-off switch	Capacity 16A /trip threshold 10A
Air inlet	In feed pressure 5.0 – 7.0 bar Flow--rate 10 l/min with Turbine 30 l/min Connection pipe diam.6x8
Suction	Suction devices powered directly from the unit: max limit of absorption, 2000 VA
Spray liquid tank capacity	1 litre
Podiatry unit top max load capacity	Kg.15
Operating temperature range	10 / 35 °C
Storage and transport temperature	-10 / 70 °C
Storage and transport humidity	10 / 90 %
Atmospheric pressure	860 / 1060 mbar
Unit weight	Kg.80 (with drawer's element)
Unit dimensions (cm)	790 x 640 x 920(H) mm
Drawer unit load capacity	Kg.3,00

The values referred to weights are indicative and vary according to the configuration requested by the customer (accessories included).

It is forbidden to use accessories other than those supplied by S.M.D.

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The symbols that appear on the equipment comply with IEC 601.1 and ISO 9687 standards.

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**GB WARNINGS****Purpose of the unit.**

The podiatry unit is intended for professional podiatrists as an aid to in-surgery operations that concern ailments of the foot and their treatment.

These units have been designed and constructed for use with patients of all ages.

Never the podiatry unit may be used in the presence of flammable anesthetic gas mixtures. Stern Podia units are not sterile devices.

The compartment with the "U.V." germicide lamp is not a sterilizing device.

When treating patients with pace-makers or similar devices the operator must take the necessary precautions: refer to the specific literature on the matter.

Do not exceed the loads indicated in the "technical characteristics" section.

Do not remove the labels on any unit/chair devices. Should these deteriorate replace them.

Do not apply any charges and devices on the tools tablet and lamp.

Is not provided using the device in emergency rooms or operating theaters.

**Power supply**

Make sure that the unit/chair electrical power supply is powered via an external differential cut-out switch with a 16A current-carrying capacity and a 10mA trip threshold.

Make sure that unit power supplies are separate.

Make sure that the podiatry unit power supply is fitted with an efficient grounding (earth) system and that wiring complies with the standards in force.

Make sure that the compressed air feed line has an external close-off valve.

Make sure that the air hoses are bled carefully and thoroughly before connecting to the podiatry unit.

**Climate**

Under extreme conditions (heat, cold, humidity) it is advisable to let a few hours pass between unpacking the equipment and using it for the first time. This precaution provides the time needed to eliminate any condensation that may have formed inside the packaging.

**Precautions**

The unit needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information contained herein.

The portable devices may affect the operation of the device.

Use only accessories that are specified in this document in order to ensure the maintenance of the EMC performance of the device.

Do not use together with other nearby devices to avoid mutual interference.

Attention: the connection of the electrical power strip leads to the creation of a SYSTEM EM, and this situation may result in a reduced level of security. About the requirements applicable to a EM SYSTEM, the ORGANIZATION RESPONSIBLE must refer to this standard.

The unit must be connected to an external air supply, whose characteristics are given in the Technical Specifications section. For the insertion of the tube (according to DIN 74324 specifies 8X1 PA11 PHL) a quick connector of 6x8 is provided.

These operations must be carried out by technicians with appropriate knowledge of the device.

For the extraction, there are no special requirements.

To preserve the condition of the buffer battery present in the unit, make sure that it does not remain disconnected from the power supply for longer than three months.

The battery replacement should be carried out by technicians with appropriate knowledge of the device.

Upon request S.M.D. provides circuit diagrams, spare parts list, descriptions, calibration instructions or other information to assist the technical support personnel during installation or repair.

Should you find any fissures or cracks on the power cord or should it be partly cut off, get SMD qualified service personnel to change it.

**S.M.D. cannot be held liable for any damages deriving from failure to observe the warnings and instructions contained in this handbook.**

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GB

## CLEANING AND DISINFECTION

To clean and disinfect the external surfaces of the podiatry unit and chair (plastic, painted, upholstered parts) use standard commercially available disinfectants and check that the listed active substances do not exceed the maximum concentrations indicated below:

- Glutaraldehyde 2%:10%
- Ethanol 96%:40%
- Formaldehyde:0.01%
- Glyoxal:0.15%
- Propanol:35%

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### WARNING!

Do not exceed the given concentrations. Do not use products containing alcohol, ammoniac, abrasive substances or benzene. Wipe disinfected areas with a wet cloth (disinfectants may attack surfaces even where diluted).

Apply products with a soft disinfectant-dampened cloth: do not spray the disinfectant directly onto equipment as it could infiltrate the unit/chair and compromise proper operation.

S.M.D. cannot be held liable for any damages deriving from failure to observe the warnings and instructions contained in this handbook.

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


GB

### SUCTION DEVICE

The suction device allows the operator to remove the organic dust (which would otherwise be inhaled) from the operating zone. This dust is produced when the micro motor is used in absence of spray.

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
#### WARNING!

Enabling of the suction device automatically disables micromotor spray.

Slot the suction tip "1" onto the hand piece/micro motor group. "2"

To start suction, remove the micro motor/suction tip set from its housing and tug the instrument.

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
**Horizontal shifting on  
"A" lever:**

-activation and adjustment  
of dynamic instruments  
rotation speed.

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
**Vertical shifting on "A"  
lever when dynamic  
instruments are no  
active:**

-activation of dynamic  
instruments chip-blower  
function.

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
**Vertical shifting on lever  
"A" when dynamic  
instruments are active:**

-spray activation both in  
continuous and sequential  
way.

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GB

### ROUTINE MAINTENANCE

To keep your Stern Podia equipment in good condition, observe the following routine maintenance schedule.

### AFTER EVERY TREATMENT SESSION OR AT LEAST DAILY

Clean and disinfect all the external surfaces of the unit.

Clean and lubricate the dynamic instruments according to the instructions provided by the instrument manufacturer.

Clean and disinfect instrument cords.

Remove the suction pipe, wash it thoroughly in warm water and neutral soap; do not sterilise it in an autoclave.

Check that everything is working properly and replace the intake filter if necessary.

### YEARLY (to be carried out by C.S.S.)

Lubricate all unit joints (movement and rotation).

Check that all safety devices are working properly.

Check that the unit is positioned in a stable manner and that the full scope movement is available.

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• Contacter le SAV.

**TROUBLESHOOTING**

Turning main switch, does not turn on the unit.

Instruments are not enabled.

No water and/or air.

Progressive reduction of fluid flowrate.

Progressive reduction of suction module/device performance.

Spray water is not fully nebulised.

Optic fibres do not function.

- Check the external differential cut-out switch.
- Check the fuses.
- Contact the CSS (Customer Support Service)
- Check the main switch.
- Contact the CSS.
- Check that pneumatic air selector is open.
- Check that external feed lines are open.
- Contact the CSS.
  
- Clean or replace the unit in feed filter cartridges (CSS)
- Check the suction filter and replace if necessary.
- Check cleaning the suction pipes and replace if necessary (CSS).
- Check that instrument spray nozzles are unobstructed; make any necessary adjustment.
- Check water level inside the tank.
- Contact the CSS
- Check that the function has been enabled.
- Contact the CSS

	<b>LISTA RICAMBI</b>	<b>PIECES DETACHEES</b>	<b>SPARE-PARTS LIST</b>
469600007	MICROMOTORE MCX	MICROMOTEUR MCX	MCX MICROMOTOR
469600006	CORDONE MICROMOTORE	CORDON MICROMOTEUR	MICROMOTOR TUBE
478990005	CORDONE TURBINA	CORDON TURBINE	TURBINE TUBE
481042000	SIRINGA COMPLETA	SERINGUE	SYRINGE
460233640	MANIGLIA CASSETTO	POIGNEE TIROIR	DRAWER HANDLE
457604011	RUOTA CON FRENO	ROUE AVEC FREIN	WHEEL WITH BRAKE

**DICHIARAZIONE DI CONFORMITA'**  
**DECLARATION DE CONFORMITE**  
**DECLARATION OF COMPLIANCE**  
**ÜBEREINSTIMMUNGSERKLÄRUNG**

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Are medical devices of Class IIa, classification rule 9 (Encl. IX of the 93/42/EEC Directive modified by the 2007/47/CE).

They follow the fundamental safety requirements as indicated in the Encl. I of the 93/42/EEC Directive (as modified by the 2007/47/CE). The quality system implemented according to the Encl. II of the 93/42/EEC Directive is certified and evaluated by the Notified Body n.0476.

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### WARNINGS ON ELECTROMAGNETIC INTERFERENCES

The device has undergone EMC (Electromagnetic Compatibility) verifications for medical devices as per EN60601-1-2.

The device needs special precautions about the EMC and must therefore be installed and commissioned following the EMC information included in this manual.

Any portable and mobile radiocommunications equipment (including also transmission peripherals such as cables and antennas) could influence the functioning of the device. Therefore, they should be placed at a distance of at least 30 cm. Failure to observe this indication could result in a malfunctioning.

Furthermore, interferences could occur in proximity of any equipment marked with the following symbol

The use of cables different from the provided ones could cause an increase of the EMISSIONS or a decrease of the IMMUNITY of the device.

The device should not be used near other electrical or electronic equipment. If necessary, verifying that the device functions correctly. In case of faults or transient noises, not only the device stops working but it also clearly indicates that. The device can be used in any building, including also residential buildings and those directly connected to the public low-voltage power supply network which supplies buildings intended for residential use.

Based on the verifications carried out, the following further warnings are provided:

The device is suited for the use in a specific electromagnetic environment. The buyer or the end-user of the device should make sure that it is used in an electromagnetic environment having the following characteristics:

Emissions test	Conformity
Radiated or conducted RF CISPR 11	Group 1 Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations and flicker emission IEC 61000-3-3	Conform

Immunity test	Test level (and conformity test) EN60601-1-2 (maximum values)
ESD EN61000-4-2	± 8kV contact discharge ± 15kV air discharge
Electrical Fast Transient (EFT) EN61000-4-4	± 2kV power line ± 1kV I/O line
Surge EN61000-4-5	± 1kV differential mode voltage ± 2kV common mode voltage
Voltage dips, short interruptions and voltage variations EN61000-4-11	0% $U_n$ during 0.5 cycle 70% $U_n$ during 25 cycles 0% $U_n$ during 5 seconds
Power-frequency magnetic field (50/60 Hz) EN61000-4-8	30A/m
Radiated RF EN61000-4-3	3V/m 80MHz ÷ 2.7GHz
Conducted RF EN61000-4-6	3V 150kHz ÷ 80MHz