I

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Bon travail !

GB

Thank you for choosing STERN PODIA products.

The equipment you have purchased has been designed and manufactured to meet your present needs and, thanks to its in-built versatility and adaptability, your future needs too.

Before using this unit, it is good practice to read this handbook carefully and keep it safely.

Doing so will provide you with a full understanding of the unit, allow you to make the best use of its potentiality and maintain top-level performance, safety, and reliability throughout its working life.

And, whatever circumstances may arise, you can count on the professionalism of your S.M.D. dealer, the expertise of his technicians and our reliability as constructors.

We wish you all the best for your work!

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GB TECHNICAL SPECIFICATIONS

Power supply	230 VAC single phase /50Hz
Podiatry unit absorption	1300 VA
Allowed mains voltage fluctuation	+/-10%
Mains fuses	n.2 (T.6,3A /250V)
Electrical protection rating	Class I -Type B
Liquids protection rating	IPX 0
Surgery electrical system	IEC 60364
Separate on-off switch	Capacity 16A /trip threshold 10A
Air line inlet	In feed pressure 5.0 – 7.0 bar Flow--rate 10 l/min with Turbine 30 l/min Connection pipe diam.6x8
Suction	Suction devices powered directly from the unit: max limit of absorption, 2000 VA
Spray liquid tank capacity	1 Lt
Tray load capacity	0,5 Kg
Lift able load	160 Kg.
Leg-rest load capacity	25 Kg
Armrest load capacity	35 Kg
Backrest load capacity	55 Kg
Operating temperature range	10 - 35 °C
Storage and transport temperature	-10°/70°C
Operating humidity range	45/95 %
Storage and transport humidity	10/90 %
Atmospheric pressure	500 - 1060 mbar
Operative group weight	Kg.140
Operative group dimensions (cm)	mm.2180 x mm.1050 x mm.1730 (H)

Weights are indicative only and may vary according to the exact configuration.

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The symbols that appear on the equipment comply with IEC 601.1 and ISO 9687 standards.

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GB WARNINGS**Purpose of the unit**

The podiatry unit is intended for professional podiatrists as an aid to in-surgery operations that concern ailments of the foot and their treatment.

These units have been designed and constructed for use with patients of all ages.

Never the podiatry unit may be used in the presence of flammable anaesthetic gas mixtures. Stern Podia units are not sterile devices.

When treating patients with pace-makers or similar devices the operator must take the necessary precautions: refer to the specific literature on the matter.

Do not exceed the loads indicated in the "technical characteristics" section.

Do not remove the labels on any unit/chair devices. Should these deteriorate replace them.

Is not provided using the device in emergency rooms or operating theaters.

Power supply

Make sure that the unit/chair electrical power supply is powered via an external differential cut-out switch with a 16A current-carrying capacity and a 10mA trip threshold.

Make sure that unit power supplies are separate.

Make sure that the podiatry unit power supply is fitted with an efficient grounding (earth) system and that wiring complies with the standards in force.

Make sure that the compressed air feed line has an external close-off valve.

Make sure that the air connections are bled carefully and thoroughly before connecting to the podiatry unit.

Climate

Under extreme conditions (heat, cold, humidity) it is advisable to let a few hours pass between unpacking the equipment and using it for the first time. This precaution provides the time needed to eliminate any condensation that may have formed inside the packaging.

Precautions

The unit needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information contained herein.

The portable devices may affect the operation of the device.

Use only accessories that are specified in this document to ensure the maintenance of the EMC performance of the device. Do not use together with other nearby devices to avoid mutual interference.

Attention to the connection of the electrical power strip leads to the creation of a SYSTEM EM, and this situation may result in a reduced level of security. About the applicable requirements SYSTEM EM the ORGANIZATION RESPONSIBLE must refer to this standard.

The unit must be connected to an external air supply, the characteristics of which are given in section Specifications: connection modes are:

inserting the plug into the outlet to be controlled by a picture corresponding to the characteristics mentioned in the technical specifications. For the insertion of the tube (according to DIN 74324 specifies 8X1 PA11 PHL), there is provided a quick connector of 6x8.

About the connection for the suction this takes place by means of a graft of diameter 30mm.

These operations must be carried out by technicians with appropriate knowledge of the device.


S.M.D. provide on request circuit diagrams, the parts list, descriptions, calibration instructions or other information to assist the technical support personnel during installation or repair.

Please note if you find any fissures or cracks on the power cord or the same was partly cut off, operate the switch activity by service personnel qualified by SMD

S.M.D. cannot be held liable for any damages deriving from failure to observe the warnings and instructions contained in this handbook.

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
GB GENERAL INFORMATION

- Keep the base constantly clean to avoid accumulation of organic deposits within the coverings and consequent stinking odours.
- Absolutely avoid pouring of liquids on the base; in case of accidental pouring, dry as soon as possible.
- Full compliance with EEC Directive 93/42 (Medical Devices) ensures maximum safety for the user, the patient, and the environment.

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❸ SAFETY

- For safety reasons, if the instrument tray arm (1) is raised accidentally a device will be triggered which will disable chair movements. To resume normal operation, bring the instrument tray arm back into the normal working position.
- EXAGONO-P units include a red emergency push button (2) which, if pressed, will bring chair movement to an immediate halt. To resume normal operation, reset via the controls.
- There is also a device which inhibits operation of the dental unit instruments if the air supply has been cut off.

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TROUBLESHOOTING

Turning main switch, does not turn on the unit.	•Check the external differential cut-out switch. •Contact the CSS (Customer Support Service)
Instruments are not enabled.	• Check the main switch. •Contact the CSS.
No water and/or air.	•Check that pneumatic air selector is open. •Check that external feed lines are open. •Contact the CSS.
Progressive reduction of fluid flow-rate.	•Clean or replace the unit in feed filter cartridges (CSS)
Progressive reduction of suction module/device performance.	•Check the suction filter and replace if necessary. •Check cleaning the suction pipes and replace if necessary (CSS).
Spray water is not fully nebulised.	•Check that instrument spray nozzles are unobstructed; make any necessary adjustment. •Check water level inside the tank. •Contact the CSS
The chair does not proceed any movement	- Check the main switch (in the back side of the chair) - Contact the CSS

	I PARTI DI RICAMBIO	F PIECES DETACHEES	GB SPARE-PARTS LIST
107020030	Bloccaggio gambiera	Blocage jambière	Leg rest blocking device
187595637	Cuscino gambiera	Sellerie jambière	Leg rest upholstery
187595875	Carter schienale	Carter dossier	Backrest cover
187595860	Cuscino seduta	Sellerie assise	Seat upholstery
187595870	Cuscino schienale	Sellerie dossier	Backrest upholstery
188778652	Appoggia testa	Têteière	Headrest
469600007	MICROMOTORE MCX	MICROMOTEUR MCX	MCX MICROMOTOR
469600006	CORDONE MICROMOTORE	CORDON MICROMOTEUR	MICROMOTOR TUBE
478990005	CORDONE TURBINA	CORDON TURBINE	TURBINE TUBE
481042000	SIRINGA COMPLETA	SERINGUE	SYRINGE

DICHIARAZIONE DI CONFORMITA'
DECLARATION DE CONFORMITE
CONFORMITY DECLARATION
KONFORMITÄTSERKLÄRUNG

S.M.D. garantisce e dichiara sotto la propria responsabilità che i dispositivi medici /
S.M.D. déclare que les dispositifs médicaux
S.M.D. declares that the medical devices / S.M.D. erklärt, dass die Medizinprodukte

Tipo / Type / Typ	Marchio/Marque/ Trademark / Brandmarke	Modello / Modèle / Model / Modell	Matricola / <u>No. de</u> fabrication /Serial No. / Seriennummer. da/de/fm/von ÷ a/a/to/bis
Riunito podologico	Stern Podia	EXAGONO	39550-00-00000 ÷ 39550-99-99999
Unité podologique			39551-00-00000 ÷ 39551-99-99999
Podiatric Unit			39552-00-00000 ÷ 39552-99-99999
Behandlungseinheit			39553-00-00000 ÷ 39553-99-99999

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Are medical devices of Class IIa, classification rule 9 (Encl. IX of the 93/42/EEC Directive modified by the 2007/47/CE). They follow the fundamental safety requirements as indicated in the Encl. I of the 93/42/EEC Directive (as modified by the 2007/47/CE). The quality system implemented according to the Encl. II of the 93/42/EEC Directive is certified and evaluated by the Notified Body n.0476.

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WARNINGS ON ELECTROMAGNETIC INTERFERENCES

The device has undergone EMC (Electromagnetic Compatibility) verifications for medical devices as per EN60601-1-2.

The device needs special precautions about the EMC and must therefore be installed and commissioned following the EMC information included in this manual.

Any portable and mobile radio communications equipment (including also transmission peripherals such as cables and antennas) could influence the functioning of the device. Therefore, they should be placed at a distance of at least 30 cm. Failure to observe this indication could result in a malfunctioning.

Furthermore, interferences could occur in proximity of any equipment marked with the following symbol

The use of cables different from the provided ones could cause an increase of the EMISSIONS or a decrease of the IMMUNITY of the device.

The device should not be used near other electrical or electronic equipment. If necessary, verifying that the device functions correctly. In case of faults or transient noises, not only the device stops working but it also clearly indicates that. The device can be used in any building, including also residential buildings and those directly connected to the public low-voltage power supply network which supplies buildings intended for residential use.

Based on the verifications carried out, the following further warnings are provided:

The device is suited for the use in a specific electromagnetic environment. The buyer or the end-user of the device should make sure that it is used in an electromagnetic environment having the following characteristics:

Emissions test	Conformity
Radiated or conducted RF	Group 1
CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations and flicker emission IEC 61000-3-3	Conform

Immunity test	Test level (and conformity test) EN60601-1-2 (maximum values)
ESD EN61000-4-2	± 8kV contact discharge ± 15kV air discharge
Electrical Fast Transient (EFT) EN61000-4-4	± 2kV power line ± 1kV I/O line
Surge EN61000-4-5	± 1kV differential mode voltage ± 2kV common mode voltage
Voltage dips, short interruptions and voltage variations EN61000-4-11	0% U _n during 0.5 cycle 70% U _n during 25 cycles 0% U _n during 5 seconds
Power-frequency magnetic field (50/60 Hz) EN61000-4-8	30A/m
Radiated RF EN61000-4-3	3V/m 80MHz ÷ 2.7GHz
Conducted RF EN61000-4-6	3V 150kHz ÷ 80MHz