



S.M.D. Srl

Via Bruno Buozzi, 41 40057 CADRIANO-GRANAROLO DELL'EMILIA (BO)  
tel. +39 051766143 fax +39 051 766043 www.smdmed.com

STERN PODIA 

P.IVA IT 02195391202 C.F./Reg.Imprese BO 02195391202 REA n.419949 Capitale sociale euro 50.000 i.v.



STERN PODIA 



STERN PODIA 

*Please note that the images present in this catalogue are for illustration and promotional purposes only. Actual product and design may slightly vary due to product enhancement, different device screen settings, the lighting in the installation location, slight differences in product finishes over time and other factors.*

*SMD Stern Podia will not accept responsibility for any color or design differences that are not factory faults.*



Steel chassis with highly resistant covers of ABS+PMMA, mounted on steel base with 4 rotary castors. Fasteners in high strength stainless steel.  
2 big lateral drawers on the left side + chest of drawers (configuration of your choice)

Worktop in Corian BTS with raised edge

Double articulated and compensated arm on three joint.

Instrument holder comprising instruments, functions and controls of electric, pneumatic and hydraulic circuits.

Pressurized container of spray water. Air filter with automatic drain of condensation.

Movable and ambidextrous Touch Screen 7" control panel with management of instruments and visualization of unit and dynamic instrument activation, as well as micromotor maximum selected speed, working speed and percentage of power air used for the turbine.

Controls for: spray activation, sequential and timer-controlled function of the sprays, optical fibers activation, rotational direction inversion (with acoustic signal) and further inlet. User set-up visualization, selection of the combined Stern Podia chair model, operating instructions and ordinary maintenance.

Wireless control of the Stern-Podia chair movements.

Safety system that prevents the simultaneous activation of the spray/suction functions.

Multifunctional "SL" foot control

Air intake: 5bar, 50 l/min.

Dimensions on floor: 50 x 78 – Operational dimensions: 156 x 103 x 155(H) cm

Weight: 115 kg

# OPERATING UNITS:



STERN PODIA

Class IIa Medical Device

Swiss and Italian  
electronic boards

Internal Transformer device  
certified as medical

Reducers to adjust both air  
and spray flow-rate


Worktop in  
Corian BTS












ABS-PMMA casings

Structure made of steel  
treated against oxidation

Wheels to guarantee  
easy movements

**Chip Blower:** air flux via direct pressure  
of the pedal lever (on SL / AS pedals)



<p style="text-align: center;"><b>DENEb</b></p>	<p style="text-align: center;"><b>LUNAR</b></p>	<p style="text-align: center;"><b>ATLAS</b></p>	<p style="text-align: center;"><b>QUASAR</b></p>
<p><b>Deneb:</b> just a few movements to change from left-handed to right-handed users and vice versa</p> 			
<p>The pedal is included in the price of the unit</p>	<p>The pedal is NOT included in the price of the unit</p> <p>Possibility to choose between:</p>	<p>The pedal is NOT included in the price of the unit</p> <p>Possibility to choose between:</p>	<p>The pedal is included in the price of the unit</p>
 <p style="text-align: center;"><b>SL Pedal</b></p>	 <p style="text-align: center;"><b>SL Pedal</b></p>  <p style="text-align: center;"><b>ON/OFF Pedal</b></p>	 <p style="text-align: center;"><b>SL Pedal</b></p>  <p style="text-align: center;"><b>ON/OFF Pedal</b></p>	 <p style="text-align: center;"><b>AS Quasar: AS Pedal</b></p>  <p style="text-align: center;"><b>A Quasar: ON/OFF Pedal</b></p>

## ON/OFF Pedal

Activation of the instruments at the pre-set speed



## SL and AS MULTIFUNCTIONAL PEDALS

allow to progressively work until the pre-set maximum speed

**SL** multifunctional pedal  
Progressive increase of the speed of  
both Micromotor and Turbine



### Horizontal shifting on "A" lever:

activation and adjustment of dynamic instruments  
(Micromotor und Turbine) rotation speed.

**AS** multifunctional pedal  
Progressive increase of the micromotor speed  
The turbine is activated to the maximum speed **only**

### Horizontal shifting on "A" lever:

Turbine activation.  
Action and adjustment of Micromotor rotation speed.




### Vertical shifting on "A" lever when dynamic instruments are not active:

activation of dynamic instruments chip-blower function.














### Vertical shifting on lever "A" when dynamic instruments are active:

spray activation both in continuous and sequential way








## SMD ADDITIONAL SECURITY SYSTEMS:

- **Hydro Stop:** air flux automatically activated any time the micromotor stops working to avoid cross contaminations
- Impossibility to work with spray and suction at the same time
- Impossibility to work if compressor is not connected to the unit
- Impossibility to work if the bottle containing water is not pressurized
- Activation of only one instrument at a time

	DENEb	LUNAR	ATLAS	QUASAR
Possibility to choose between			F	F
<b>B: spring tension instruments</b> 	B / F	B / F		
<b>F: hanging hoses instruments</b> 				
<b>Control Panel</b>	7" Touch Screen 	5" Touch Screen 	LCD 	LED light 
<b>Wireless System</b> for controlling the chair movements				-
<b>Maximum number of Instruments</b> (dynamic instruments are either Micromotor or Turbine)	<b>4</b> (3 of them can be <b>dynamic</b> instruments)	<b>4</b> (3 of them can be <b>dynamic</b> instruments)	<b>4</b> (3 of them can be <b>dynamic</b> instruments)	<b>4</b> (3 of them can be <b>dynamic</b> instruments)



	DENEb	LUNAR	ATLAS	QUASAR
Worktop in <b>Corian BTS</b>	X	X	X	-
Possibility to choose <b>drawers configuration</b> 	X	X	X	-
Intermittent or continuous <b>Spray function</b> → Available in <b>SMD operating unit only</b> 	X	X	X	-
<b>ROTATION DIRECTION PUSH-BUTTONS</b>  	X	X	X	X
<b>OPTICAL FIBRES SWITCH ON/OFF PUSH-BUTTONS</b> 	X	X	X	-





### WITHOUT HANDLES

Drawers mounted on "push" slides that allow the opening of the drawer with a slight pressure on the front panel.



### WHITE HANDLES

Drawers mounted on soft closing slides and equipped with a newly designed **aluminum handle painted in white**.



### COLORED HANDLES

Drawers mounted on soft closing slides and equipped with a new design aluminum handle **painted in a color of your choice from Stern Podia selection**.

## 3 NEW TYPES OF DRAWERS

# Instruments at work ...



## MICROMOTORS

indicated for all abrasion treatments on nails laminas and skin thickness

Application for Electric micro motor BIEN AIR MX-2 Brush-less, 40Krpm, with Optical Fibres LED Technology, integrated spray, hidrostop device, chip-blower (HAND PIECE NOT INCLUDED).



Application for Electric micro motor BIEN AIR MX- i Brush-less 40Krpm, self-ventilated (HAND PIECE NOT INCLUDED).



Application for Electric micro motor BIEN AIR MC40, 40Krpm, integrated spray, hidrostop device, chip-blower (HAND PIECE NOT INCLUDED).



Pneumatic micromotor Bien-Air Aquilon with Midwest coupling, integrated spray (HAND PIECE NOT INCLUDED).



## HAND-PIECES

STERN PODIA Hand-piece QUARK, optical fibers with LED technology, integrated spray






























Hand-piece Bien Air PM1:1 Micro Series with integrated spray













STERN PODIA Hand-piece QUARK, integrated spray



	Good
	Not suggested
	Not indicated

MICROMOTORS	Type	Spray	Thermic protection	Optical fibres	Torque at 20.000 rpm	Torque Compensation	Small Burrs (-5mm)			Big Burrs (+ 5mm)		
							Working time			Working time		
							2 h/d	2-4 h/d	+ 4 h/d	-2 h/d	2-4 h/d	+ 4 h/d
<b>Bien Air MX-2</b>	Brushless	YES	YES I	YES	3,5 Ncm	YES						
<b>Bien Air MX-i Self-ventilated</b>	Brushless	-	YES	-	5 Ncm	YES						
<b>Bien Air MC40</b>	With Brushes	YES	-	-	2,3 Ncm	-						
<b>Bien Air Aquilon</b>	Pneumatic	YES	-	-	2,5 Ncm	-						

HANDPIECES	Bien Air MX2	Bien Air MX-i	Bien Air MC40	Bien Air Aquilon
<b>Bien Air PM 1:1 Micro Series</b>				
<b>Stern Podia Quark LED technology</b>				
<b>Stern Podia Quark</b>				

“Straight” 3-way syringe with air, water, and spray



Integrated suction duct with stop device for spring tension arm and coaxial nozzle to be applied to the micro motor. Suction effect prolonged and timer-controlled. - indicated to work with dry or mixed technique, suitable for all micromotors in the Price List





## TURBINE

indicated for all abrasion treatments on nails corners and skin thickness.

Arrangement for turbine with Midwest coupling, optical fibres lightening system, spray, chip-blower, hidrostop device (hand piece not included).

Arrangement necessary for the functioning of the turbine

Turbine TKD MIKRA R model, Midwest connection for rounded adaptor



Rounded Adaptor Midwest connection



Turbine TKD BRAVIA LB model, with LED light, Midwest connection for rounded adaptor



Rounded Adaptor Midwest connection



### ASPI-COMP Suction motor

The suction motor is housed in a proper sound-proof steel cabinet, with silencer, switch and protection fuses.

Suction filter included  
Sucked air: 60 m<sup>3</sup>/h.



### ASPI-COMP Suction motor with suction arm



### ASPI-COMP integrated in the operating unit

(not available for Quasar and Exagono models)



# Operating units and group accessories



### **Genesis 24 air compressor**

Oil-less 2-cylinders motor with soundproofing case membrane dryer, air-condensation filter 0,01micron

- Power 1,1KW
- Air delivery 120 l/min
- Tank 24 liter
- Noise 50 dB(A)
- Dimensions 530x500x800
- Weight 67 Kg



### **Post for operating lights**

Available in three different heights:

- 530 mm
- 340 mm
- 220 mm



### **SUNLED led lamp**

7 led x 3 Watts (used 70%)  
Color temperature 6.300° K  
Power supply 24 Volts  
15 watts total  
Double articulated and compensated arm



# Operating units and group accessories


**Tray for  
instrument  
holder**



**Inserts for drawers**



**Additional chest of drawer**



Operating units and group accessories



Operating units and group accessories

STERN PODIA **S.M.D. Srl**

Via Bruno Buozzi, 41 40057 CADRIANO-GRANAROLO DELL'EMILIA (BO)  
tel. +39 051766143 fax +39 051 766043 [www.smdmed.com](http://www.smdmed.com)

STERN PODIA 

P.IVA IT 02195391202 C.F./Reg.Imprese BO 02195391202 REA n.419949 Capitale sociale euro 50.000 i.v.