



STERN PODIA 

	<p>S.M.D. Srl Via Bruno Buozzi, 41 40057 CADRIANO-GRANAROLO DELL'EMILIA (BO) tel. +39 051766143 fax +39 051 766043 www.smdmed.com</p>	STERN PODIA 
<p>P.IVA IT 02195391202 C.F./Reg.Imprese BO 02195391202 REA n.419949 Capitale sociale euro 50.000 i.v.</p>		

*Le immagini presenti in questo catalogo sono inserite a scopo illustrativo ed hanno il solo scopo di presentare il prodotto per la vendita, senza alcuna garanzia o impegno, da parte di SMD Stern Podia, circa l'esatta corrispondenza dell'immagine raffigurata con il prodotto reale
SMD Stern Podia non accetterà responsabilità per differenze in termini di colore o design che non siano dovute ad errori di produzione.*

CAB X2

- **Hermetic closure:** no dust exit - no dangerous inhalations
- **No vacuum:** energy saving and no noise
- **Great visibility:** wide high-strength tempered glass and 24 LEDs light
- **Metallic structure:** low electrostatic properties - no dust attraction
- **Specific Height** chosen for operators of 1,55 m y 1,90 m height.






- **Big working area**
- **Easiness of work:** a pedal activates the abrasive wheels, no need of safety systems because the operator presence is needed to activate the machine



CAB X2

- **Innovative system with two independent motors:** working field always in a central position. A switch on the front panel allows selecting the desired motor to alternatively activate the abrasive wheel or the finishing brush.
- **Single shaft between each motor and the relevant abrasive wheel:** no vibrations during work and no noise for a very precise work.
- **Sleeves that protects operator and avoid the exit of dust:** hermetic closure



STERN PODIA **S.M.D. Srl**

Via Bruno Buozzi, 41 40057 CADRIANO-GRANAROLO DELL'EMILIA (BO)
tel. +39 051766143 fax +39 051 766043 www.smdmed.com

STERN PODIA 

P.IVA IT 02195391202 C.F./Reg.Imprese BO 02195391202 REA n.419949 Capitale sociale euro 50.000 i.v.